



# Thermo Scientific™ Sensititre™ System

For identification & antimicrobial susceptibility testing

# Facilitate optimal patient outcomes and respond more effectively to emerging resistance

## Thermo Scientific™ Sensititre™ System

With a proven history of AST accuracy, the Sensititre™ System delivers definitive MIC results for both commonly prescribed and novel antimicrobials.

While other automated AST systems may rely on generating MIC values based on extrapolated data calculated from historical data sets, the Sensititre™ System utilizes gold standard-level MIC accuracy, providing superior quality and reproducibility for definitive results the first time.

### More antimicrobials, more testing options

The Sensititre™ System combines manual, semi-automated, and fully automated testing on a single AST platform using broth microdilution (BMD) - the gold standard for AST accuracy. With access to over 300 antimicrobials, including more than 40 veterinary-specific antimicrobials, as well as offering the latest-to-market on standard and custom MIC plates,

the Sensititre™ System allows you to consolidate confirmatory tests and can offer plates that meet US FDA, CLSI, and/or EUCAST breakpoint requirements. Moreover, the Sensititre™ System provides a standardized AMR surveillance tool to support public health and national reference laboratories, enabling you to comply with government surveillance mandates while streamlining workflow. Harmonize your surveillance AST with the method of choice for global AMR programs, including:

- U.S. National Antimicrobial Resistance Monitoring System (NARMS), coordinated via FDA/CVM, USDA, and the CDC
- U.S. CDC Antibiotic Resistance Lab Network
- EU Monitoring System of Zoonoses

## Powerful, comprehensive solutions for common and complex pathogens

Antimicrobial resistance (AMR) has reached global health crisis levels. This, combined with a lack of comprehensive AMR data to react to AMR trends, poses challenges to many microbiology labs, especially post-COVID-19. Robust antimicrobial susceptibility testing (AST) is critical in the fight against AMR. It aids in optimizing treatment for critically ill patients, supporting enhanced patient outcomes, stewardship initiatives, and worldwide AMR surveillance in human, animal, and environmental sectors.

It is also instrumental in safeguarding our existing antimicrobials against escalating resistance.

### Empowering your expertise: from test to results

From manual disc diffusion to automated instrumentation with definitive minimum inhibitory concentration (MIC) results, our Thermo Scientific AST portfolio can be used for a wide array of applications, from common, everyday AST to specialty applications such as yeast and mycobacteria, with gold standard-level results<sup>1</sup>. Our AST solutions can help clinicians, veterinarians, food scientists, and public health officials determine pathogen antimicrobial susceptibility patterns, while supporting global and regional antimicrobial stewardship.

# Consolidate all AST on a single platform

Meeting the unique testing demands of all laboratory sizes and volumes, choose from a variety of manual and automated options

## Select Plate

### Sensititre™ Standard AST Plates:

Choose from our wide selection of standard plates to suit your application, including Gram-positive, Gram-negative, fastidious, mycobacteria, and yeast formats, as well as veterinary and surveillance testing formats.



### Sensititre™ Custom AST Plates:

Design your own plate from our selection of over 300 antimicrobials and over 40 veterinary-specific antimicrobials. We offer one of the widest, most up-to-date selections of antimicrobials, available in wide dilution ranges.



## Inoculum

### Thermo Scientific™ Sensititre™ Nephelometer:

A simple solution for inoculum density measurements and standardized inoculation preparation.



### Thermo Scientific™ Sensititre™ 8-Channel Programmable Pipette:

Quickly and accurately inoculate microtitre plates, with enhanced ergonomics.



### Thermo Scientific™ Sensititre™ AIM™ Automated Inoculation Delivery System:

Automatically doses Sensititre™ plates and eliminates skipped wells and costly repeat tests.



## Incubate

### Thermo Scientific™ Sensititre™ ARIS HiQ™ System:

Individually incubates 100 MIC, breakpoint, or identification plates, ensuring optimal growth conditions and eliminating offline or supplementary testing.



### Thermo Scientific™ Sensititre™ OptiRead™ Automated Fluorometric Plate Reading System:

Utilizes fluorescence detection technology to automate Sensititre™ plate reading, delivering fast, accurate results, and directly linking to the Sensititre™ SWIN Software System to automate interpretation and result reporting.



### Thermo Scientific™ Sensititre™ Vizion™ Digital MIC Viewing System:

Captures and stores easy-to-read digital plate images for optimized manual reading and traceability and connects to the Sensititre™ SWIN Software System for automated interpretation and reporting of results.



### Thermo Scientific™ Sensititre™ Manual Viewer:

Perform simple visual reads of your 96-well microtitre plates with our mirrored viewbox.



## Interpret

### Thermo Scientific™ Sensititre™ SWIN™ Software System:

Consolidates Sensititre™ results from manual and automated reading options on a single software platform.




### Sensititre™ SWIN Epidemiology Module:

Provides comprehensive reporting of AST results to help detect, monitor, and track local antibiotic resistance patterns to facilitate decision-making and support antimicrobial stewardship programs.



# Sensititre™ GNID Plate




Put 3-5 colonies into H<sub>2</sub>O to measure a 0.5 McFarland using the Nephelometer

Inoculate plate with 50 µL volume per well of the suspension using the Sensititre™ AIM or Multi-Channel Pipette

Add mineral oil to wells 1A, 2A, 5A, 6A, 9A and 10A

Seal Sensititre™ plate and incubate at 34-36°C in a non-CO<sub>2</sub> incubator or Sensititre™ ARIS HiQ for 5-18 hours


Read automatically with ARIS HiQ or OptiRead



## Antimicrobics

|             |                                                      |
|-------------|------------------------------------------------------|
| <b>FR1</b>  | Lysine AMC                                           |
| <b>UREA</b> | Urea                                                 |
| <b>ORN</b>  | Ornithine                                            |
| <b>SORB</b> | SORB Sorbitol                                        |
| <b>XYL</b>  | Xylose                                               |
| <b>FR12</b> | 2-Acetamido-2-deoxyglucoside MU + α-L-Arabinoside MU |
| <b>SUC</b>  | Sucrose                                              |
| <b>FR9</b>  | β-D-glucoside MU                                     |
| <b>FR3</b>  | Phosphate MU                                         |
| <b>TREH</b> | Trehalose                                            |
| <b>FR8</b>  | bis-Phosphate MU                                     |
| <b>MANN</b> | Mannitol                                             |
| <b>MALT</b> | Maltose                                              |
| <b>FR4</b>  | α-D-galactoside MU                                   |
| <b>INOS</b> | Inositol                                             |
| <b>FR10</b> | β-D-glucoside MU                                     |
| <b>FR5</b>  | Proline AMC                                          |
| <b>FRUC</b> | Fructose                                             |
| <b>ESCU</b> | Esculin                                              |
| <b>ATOL</b> | Arabitol                                             |
| <b>ANOS</b> | Arabinose                                            |
| <b>LYS</b>  | Lysine                                               |
| <b>TDA</b>  | TDA                                                  |
| <b>RAFF</b> | Raffinose                                            |
| <b>FR7</b>  | γ-Glutamine MU                                       |
| <b>ARG</b>  | Arginine                                             |
| <b>FR6</b>  | α-D-galactoside MU                                   |
| <b>CELL</b> | Cellobiose                                           |
| <b>MALN</b> | Malonate                                             |
| <b>PYR</b>  | Pyruvate                                             |
| <b>CIT</b>  | Citrate                                              |

# Sensititre™ GPID Plate



Put 3-5 colonies into H<sub>2</sub>O to measure a 0.5 McFarland using the Nephelometer

Inoculate plate with 50 µL volume per well of the suspension using the Sensititre™ AIM or Multi-Channel Pipette

Add mineral oil to wells 1A, 5A, and 9A

Seal Sensititre™ plate and incubate at 34-36°C in a non-CO<sub>2</sub> incubator or Sensititre™ ARIS HiQ for 18-24 hours

Read automatically with ARIS HiQ or OptiRead

|   | 1     | 2     | 3     | 4     | 5     | 6     | 7     | 8     | 9     | 10    | 11    | 12    |
|---|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| A | UREA  | ESCU  | ARG   | FR 13 | UREA  | ESCU  | ARG   | FR 13 | UREA  | ESCU  | ARG   | FR 13 |
| B | FR 16 | RHAM  | FR 15 | MANN  | FR 16 | RHAM  | FR 15 | MANN  | FR 16 | RHAM  | FR 15 | MANN  |
| C | FR 14 | TREH  | FR 17 | MALT  | FR 14 | TREH  | FR 17 | MALT  | FR 14 | TREH  | FR 17 | MALT  |
| D | FR 18 | FR 19 | FR 20 | FR 21 | FR 18 | FR 19 | FR 20 | FR 21 | FR 18 | FR 19 | FR 20 | FR 21 |
| E | FR 22 | FR 23 | GLYC  | FR 24 | FR 22 | FR 23 | GLYC  | FR 24 | FR 22 | FR 23 | GLYC  | FR 24 |
| F | FR 25 | GLUC  | SUC   | FR 26 | FR 25 | GLUC  | SUC   | FR 26 | FR 25 | GLUC  | SUC   | FR 26 |
| G | FR 27 | B-MG  | FR 28 | FR 29 | FR 27 | B-MG  | FR 28 | FR 29 | FR 27 | B-MG  | FR 28 | FR 29 |
| H | SORB  | FR 30 | FR 31 | FR 32 | SORB  | FR 30 | FR 31 | FR 32 | SORB  | FR 30 | FR 31 | FR 32 |

## Antimicrobics

|             |                             |
|-------------|-----------------------------|
| <b>UREA</b> | Urea                        |
| <b>ESCU</b> | Esculin                     |
| <b>ARG</b>  | Arginine                    |
| <b>FR13</b> | 4MUβ- D - Galactopyranoside |
| <b>FR16</b> | 4MUβ- D - Ribofuranoside    |
| <b>RHAM</b> | Rhamnose                    |
| <b>FR15</b> | Alanine 7AMC                |
| <b>MANN</b> | Mannitol                    |
| <b>FR14</b> | 4MUβ- D - Glucopyranoside   |
| <b>TREH</b> | Trehalose                   |
| <b>FR17</b> | D - Alanine 7AMC            |
| <b>MALT</b> | Maltose                     |
| <b>FR18</b> | 4MUβ- D - Mannopyranoside   |
| <b>FR19</b> | Ornithine 7AMC              |
| <b>FR20</b> | Arginine 7AMC               |
| <b>FR21</b> | 4MUβ- D - Glucuronide       |
| <b>FR22</b> | 4MU α - D - Glucopyranoside |
| <b>FR23</b> | Cysteine 7AMC               |
| <b>GLYC</b> | Glycerol                    |
| <b>FR24</b> | Threonine 7AMC              |
| <b>FR25</b> | Methionine 7AMC             |
| <b>GLUC</b> | Glucoside                   |
| <b>SUC</b>  | Sucrose                     |
| <b>FR26</b> | Proline 7AMC                |
| <b>FR27</b> | Serine 7AMC                 |
| <b>B-MG</b> | β-Methyl Glucoside          |
| <b>FR28</b> | Citrulline 7AMC             |
| <b>FR29</b> | Pyroglutamate 7AMC          |
| <b>SORB</b> | Sorbitol                    |
| <b>FR30</b> | Tyrosine 7AMC               |
| <b>FR31</b> | Leucine 7AMC                |
| <b>FR32</b> | Valine 7AMC                 |

# Sensititre™ Gram Positive GPALL1F Plate

| Intended use                         | Organism                              | Use and methodology | Inoculum preparation                                                | Broth Type                               | Read method                                                                                      | Pack Size   |
|--------------------------------------|---------------------------------------|---------------------|---------------------------------------------------------------------|------------------------------------------|--------------------------------------------------------------------------------------------------|-------------|
| Antimicrobial Susceptibility Testing | Non-fastidious Gram positive isolates | IVD-FDA (CLSI)      | Sensititre™ Sterile Water (T3339)<br>0.5 McFarland Standard (E1041) | Sensititre™ Mueller Hinton Broth (T3462) | Sensititre™ ARIS HIQ<br>Sensititre™ OptiRead<br>Sensititre™ Vizion<br>Sensititre™ Manual Viewbox | 10 Isolates |

Put 3-5 colonies into H<sub>2</sub>O to measure a 0.5 McFarland using the Nephelometer, mix 10 µL or 30 µL\* of suspension into MHB

Inoculate plate with 50 µL volume per well of the suspension using the Sensititre™ AIM or Multi-Channel Pipette

Seal Sensititre™ plate and incubate at 34-36°C in a non-CO<sub>2</sub> incubator or Sensititre™ ARIS HIQ for 18-24 hours

Read automatically with ARIS HIQ or OptiRead, semi-automatically with Vizion or manually with Manual Viewbox


|   | 1           | 2        | 3        | 4        | 5         | 6        | 7        | 8        | 9       | 10       | 11     | 12     |
|---|-------------|----------|----------|----------|-----------|----------|----------|----------|---------|----------|--------|--------|
| A | CHL 2       | CHL 4    | CHL 8    | CHL 16   | ERY 0.25  | ERY 0.5  | ERY 1    | ERY 2    | ERY 4   | CLI 0.5  | CLI 1  | CLI 2  |
| B | DAP 0.5     | DAP 1    | DAP 2    | DAP 4    | OXA+ 0.25 | OXA+ 0.5 | OXA+ 1   | OXA+ 2   | OXA+ 4  | STR 1000 | DT1    | DT2    |
| C | GEN 2       | GEN 4    | GEN 8    | GEN 16   | AMP 0.12  | AMP 0.25 | AMP 0.5  | AMP 1    | AMP 2   | AMP 4    | AMP 8  | FOXS 6 |
| D | LZD 1       | LZD 2    | LZD 4    | LZD 8    | PEN 0.06  | PEN 0.12 | PEN 0.25 | PEN 0.5  | PEN 1   | PEN 2    | PEN 4  | PEN 8  |
| E | RIF 0.5     | RIF 1    | RIF 2    | RIF 4    | VAN 0.25  | VAN 0.5  | VAN 1    | VAN 2    | VAN 4   | VAN 8    | VAN 16 | VAN 32 |
| F | SXT 0.5/9.5 | SXT 1/19 | SXT 2/38 | SXT 4/76 | LEVO 0.25 | LEVO 0.5 | LEVO 1   | LEVO 2   | LEVO 4  | CIP 1    | CIP 2  | POS    |
| G | SYN 0.5     | SYN 1    | SYN 2    | SYN 4    | TGC 0.03  | TGC 0.06 | TGC 0.12 | TGC 0.25 | TGC 0.5 | NIT 32   | NIT 64 | POS    |
| H | TET 2       | TET 4    | TET 8    | TET 16   | MXF 0.25  | MXF 0.5  | MXF 1    | MXF 2    | MXF 4   | GEN 500  | NEG    | POS    |

## Antimicrobics

|             |                               |
|-------------|-------------------------------|
| <b>AMP</b>  | Ampicillin                    |
| <b>FOXS</b> | Cefoxitin screen              |
| <b>CHL</b>  | Chloramphenicol               |
| <b>CIP</b>  | Ciprofloxacin                 |
| <b>CLI</b>  | Clindamycin                   |
| <b>DT1</b>  | D Test 1                      |
| <b>DT2</b>  | D Test 2                      |
| <b>DAP</b>  | Daptomycin                    |
| <b>ERY</b>  | Erythromycin                  |
| <b>GEN</b>  | Gentamicin                    |
| <b>LEVO</b> | Levofloxacin                  |
| <b>LZD</b>  | Linezolid                     |
| <b>MXF</b>  | Moxifloxacin                  |
| <b>NEG</b>  | Negative control              |
| <b>NIT</b>  | Nitrofurantoin                |
| <b>OXA+</b> | Oxacillin + 2% NaCl           |
| <b>PEN</b>  | Penicillin                    |
| <b>SYN</b>  | Quinupristin/Dalfopristin     |
| <b>RIF</b>  | Rifampin                      |
| <b>STR</b>  | Streptomycin                  |
| <b>TET</b>  | Tetracycline                  |
| <b>TGC</b>  | Tigecycline                   |
| <b>SXT</b>  | Trimethoprim/Sulfamethoxazole |
| <b>VAN</b>  | Vancomycin                    |
| <b>POS</b>  | Positive control              |

\*For aid in detection of resistance mechanisms for Gram positive organisms

# Sensititre™ Vancomycin KRVANCO Plate




## Antimicrobics

- NEG Negative Control
- VAN Vancomycin
- POS Positive Control

# Sensititre™ Gram Negative DKMGN Plate

with Colistin, Ceftazidime/Avibactam and Ceftolozane/Tazobactam



Put 3-5 colonies into H<sub>2</sub>O to measure a 0.5 McFarland using the Nephelometer, mix 1 µl\*, 10 µl\*\*, or 30 µl\*\*\* of suspension into MHb

Inoculate plate with 50 µl volume per well of the suspension using the Sensititre™ AIM or Multi-Channel Pipette

Seal Sensititre™ plate and incubate at 34-36°C in a non-CO<sub>2</sub> incubator or Sensititre™ ARIS HiQ for 18-24 hours

Read semi-automatically with Vizion or manually with Manual Viewbox


|   | 1            | 2            | 3            | 4            | 5            | 6            | 7            | 8           | 9            | 10          | 11           | 12           |
|---|--------------|--------------|--------------|--------------|--------------|--------------|--------------|-------------|--------------|-------------|--------------|--------------|
| A | MERO<br>0.12 | MERO<br>0.25 | MERO<br>0.5  | MERO<br>1    | MERO<br>2    | MERO<br>4    | MERO<br>8    | MERO<br>16  | AMI<br>4     | AMI<br>8    | AMI<br>16    | AMI<br>32    |
| B | GEN<br>0.5   | GEN<br>1     | GEN<br>2     | GEN<br>4     | GEN<br>8     | AZT<br>0.5   | AZT<br>1     | AZT<br>2    | AZT<br>4     | AZT<br>8    | AZT<br>16    | AZT<br>32    |
| C | CIP<br>0.06  | CIP<br>0.12  | CIP<br>0.25  | CIP<br>0.5   | CIP<br>1     | CIP<br>2     | P/T4<br>1/4  | P/T4<br>2/4 | P/T4<br>4/4  | P/T4<br>8/4 | P/T4<br>16/4 | P/T4<br>32/4 |
| D | AUGC<br>4/2  | AUGC<br>8/2  | AUGC<br>16/2 | AUGC<br>32/2 | AUGC<br>64/2 | C/T<br>0.5/4 | C/T<br>1/4   | C/T<br>2/4  | C/T<br>4/4   | C/T<br>8/4  | C/T<br>16/4  | C/T<br>32/4  |
| E | COL<br>0.25  | COL<br>0.5   | COL<br>1     | COL<br>2     | COL<br>4     | COL<br>8     | FOT<br>0.5   | FOT<br>1    | FOT<br>2     | FOT<br>4    | FOT<br>8     | TOB<br>1     |
| F | TGC<br>0.25  | TGC<br>0.5   | TGC<br>1     | TGC<br>2     | TGC<br>4     | SXT<br>1/19  | SXT<br>2/38  | SXT<br>4/76 | SXT<br>8/152 | TOB<br>2    | TOB<br>4     | TOB<br>8     |
| G | TAZ<br>0.5   | TAZ<br>1     | TAZ<br>2     | TAZ<br>4     | TAZ<br>8     | TAZ<br>16    | CZA<br>0.5/4 | CZA<br>1/4  | CZA<br>2/4   | CZA<br>4/4  | CZA<br>8/4   | CZA<br>16/4  |
| H | IMI<br>0.5   | IMI<br>1     | IMI<br>2     | IMI<br>4     | IMI<br>8     | IMI<br>16    | ETP<br>0.12  | ETP<br>0.25 | ETP<br>0.5   | ETP<br>1    | ETP<br>2     | POS          |

## Antimicrobics

|             |                                          |
|-------------|------------------------------------------|
| <b>MERO</b> | Meropenem                                |
| <b>GEN</b>  | Gentamicin                               |
| <b>CIP</b>  | Ciprofloxacin                            |
| <b>AUGC</b> | Amoxicillin / Clavulanic acid constant 2 |
| <b>COL</b>  | Colistin                                 |
| <b>TGC</b>  | Tigecycline                              |
| <b>TAZ</b>  | Ceftazidime                              |
| <b>IMI</b>  | Imipenem                                 |
| <b>AZT</b>  | Aztreonam                                |
| <b>C/T</b>  | Ceftolozane / Tazobactam 4               |
| <b>SXT</b>  | Trimethoprim / Sulfamethoxazole          |
| <b>P/T4</b> | Piperacillin / Tazobactam constant 4     |
| <b>FOT</b>  | Cefotaxime                               |
| <b>CZA</b>  | Ceftazidime / Avibactam                  |
| <b>ETP</b>  | Ertapenem                                |
| <b>AMI</b>  | Amikacin                                 |
| <b>TOB</b>  | Tobramycin                               |
| <b>POS</b>  | Positive control                         |

\*For Proteus spp. \*\*For Enterobacteriaceae and Non-Enterobacteriaceae. \*\*\*For aid in detection of resistance mechanisms for Enterobacteriaceae and Non-Enterobacteriaceae.

# Sensititre™ Gram Negative GN4F Plate



Put 3-5 colonies into H<sub>2</sub>O to measure a 0.5 McFarland using the Nephelometer, mix 1 µl, 10 µl\*\*, or 30 µl\*\*\* of suspension into MHB

Inoculate plate with 50 µl volume per well of the suspension using the Sensititre™ AIM or Multi-Channel Pipette

Seal Sensititre™ plate and incubate at 34-36°C in a non-CO<sub>2</sub> incubator or Sensititre™ ARIS HIQ for 18-24 hours

Read automatically with ARIS HIQ or OptiRead, semi-automatically with Vizion or manually with Manual Viewbox




|   | 1             | 2            | 3         | 4          | 5           | 6          | 7           | 8         | 9            | 10        | 11         | 12         |
|---|---------------|--------------|-----------|------------|-------------|------------|-------------|-----------|--------------|-----------|------------|------------|
| A | AMI<br>8      | TGC<br>1     | LEVO<br>1 | DOR<br>0.5 | ETP<br>0.25 | IMI<br>0.5 | MERO<br>0.5 | FAZ<br>1  | TAZ<br>1     | AZT<br>1  | FEP<br>4   | AXO<br>0.5 |
| B | AMI<br>16     | TGC<br>2     | LEVO<br>2 | DOR<br>1   | ETP<br>0.5  | IMI<br>1   | MERO<br>1   | FAZ<br>2  | TAZ<br>2     | AZT<br>2  | FEP<br>8   | AXO<br>1   |
| C | AMI<br>32     | TGC<br>4     | LEVO<br>4 | DOR<br>2   | ETP<br>1    | IMI<br>2   | MERO<br>2   | FAZ<br>4  | TAZ<br>4     | AZT<br>4  | FEP<br>16  | AXO<br>2   |
| D | P/T4<br>8/4   | TGC<br>8     | LEVO<br>8 | DOR<br>4   | ETP<br>2    | IMI<br>4   | MERO<br>4   | FAZ<br>8  | TAZ<br>8     | AZT<br>8  | FEP<br>32  | AXO<br>4   |
| E | P/T4<br>16/4  | TIM2<br>8/2  | NIT<br>32 | MIN<br>1   | ETP<br>4    | IMI<br>8   | MERO<br>8   | FAZ<br>16 | TAZ<br>16    | AZT<br>16 | CIP<br>0.5 | AXO<br>8   |
| F | P/T4<br>32/4  | TIM2<br>16/2 | NIT<br>64 | MIN<br>2   | ETP<br>8    | PIP<br>16  | GEN<br>2    | TOB<br>2  | A/S2<br>4/2  | AMP<br>8  | CIP<br>1   | AXO<br>16  |
| G | P/T4<br>64/4  | TIM2<br>32/2 | TET<br>4  | MIN<br>4   | SXT<br>2/38 | PIP<br>32  | GEN<br>4    | TOB<br>4  | A/S2<br>8/4  | AMP<br>16 | CIP<br>2   | AXO<br>32  |
| H | P/T4<br>128/4 | TIM2<br>64/2 | TET<br>8  | MIN<br>8   | SXT<br>4/76 | PIP<br>64  | GEN<br>8    | TOB<br>8  | A/S2<br>16/8 | POS       | POS        | POS        |

## Antimicrobics

|      |                                          |
|------|------------------------------------------|
| AMI  | Amikacin                                 |
| P/T4 | Piperacillin / Tazobactam constant 4     |
| TGC  | Tigecycline                              |
| TIM2 | Ticarcillin / Clavulanic acid constant 2 |
| LEVO | Levofloxacin                             |
| NIT  | Nitrofurantoin                           |
| TET  | Tetracycline                             |
| DOR  | Doripenem                                |
| MIN  | Minocycline                              |
| ETP  | Ertapenem                                |
| SXT  | Trimethoprim / Sulfamethoxazole          |
| IMI  | Imipenem                                 |
| PIP  | Piperacillin                             |
| MERO | Meropenem                                |
| GEN  | Gentamicin                               |
| FAZ  | Ceftazidime                              |
| TOB  | Tobramycin                               |
| TAZ  | Ceftazidime                              |
| A/S2 | Ampicillin / Sulbactam 2:1 ratio         |
| AZT  | Aztreonam                                |
| AMP  | Amicillin                                |
| FEP  | Cefepime                                 |
| CIP  | Ciprofloxacin                            |
| AXO  | Ceftriaxone                              |
| POS  | Positive control                         |

\*For Proteus spp. \*\*For Enterobacteriaceae and Non-Enterobacteriaceae. \*\*\*For aid in detection of resistance mechanisms for Enterobacteriaceae and Non-Enterobacteriaceae.

# Sensititre™ Gram Negative FRCOL Plate



## Antimicrobics

COL Colistin  
POS Positive Control

\*For Proteus spp. \*\*For Enterobacteriaceae and Non-Enterobacteriaceae. \*\*\*For aid in detection of resistance mechanisms for Enterobacteriaceae and Non-Enterobacteriaceae.

# Sensititre™ Gram Negative EURGNCOL Plate

with Colistin, Ceftazidime/Avibactam and Ceftolozane/Tazobactam

| Intended use                         | Organism                              | Use and methodology | Inoculum preparation                                                | Broth Type                               | Read method                                      | Pack Size   |
|--------------------------------------|---------------------------------------|---------------------|---------------------------------------------------------------------|------------------------------------------|--------------------------------------------------|-------------|
| Antimicrobial Susceptibility Testing | Non-fastidious Gram Negative isolates | IVD (EUCAST)        | Sensititre™ Sterile Water (T3339)<br>0.5 McFarland Standard (E1041) | Sensititre™ Mueller Hinton Broth (T3462) | Sensititre™ Vizion<br>Sensititre™ Manual Viewbox | 30 Isolates |

Put 3-5 colonies into H<sub>2</sub>O to measure a 0.5 McFarland using the Nephelometer, mix 1 µl\*, 10 µl\*\*, or 30 µl\*\*\* of suspension into MHb

Inoculate plate with 50 µl volume per well of the suspension using the Sensititre™ AIM or Multi-Channel Pipette

Seal Sensititre™ plate and incubate at 34-36°C in a non-CO<sub>2</sub> incubator or Sensititre™ ARIS HiQ for 18-24 hours

Read semi-automatically with Vizion or manually with Manual Viewbox

|   | 1           | 2             | 3            | 4            | 5           | 6             | 7            | 8            | 9           | 10            | 11           | 12           |
|---|-------------|---------------|--------------|--------------|-------------|---------------|--------------|--------------|-------------|---------------|--------------|--------------|
| A | COL<br>0.25 | P/T4<br>4/4   | C/T<br>4/4   | POS          | COL<br>0.25 | P/T4<br>4/4   | C/T<br>4/4   | POS          | COL<br>0.25 | P/T4<br>4/4   | C/T<br>4/4   | POS          |
| B | COL<br>0.5  | P/T4<br>8/4   | C/T<br>8/4   | MERO<br>0.25 | COL<br>0.5  | P/T4<br>8/4   | C/T<br>8/4   | MERO<br>0.25 | COL<br>0.5  | P/T4<br>8/4   | C/T<br>8/4   | MERO<br>0.25 |
| C | COL<br>1    | P/T4<br>16/4  | CZA<br>1/4   | MERO<br>0.5  | COL<br>1    | P/T4<br>16/4  | CZA<br>1/4   | MERO<br>0.5  | COL<br>1    | P/T4<br>16/4  | CZA<br>1/4   | MERO<br>0.5  |
| D | COL<br>2    | P/T4<br>32/4  | CZA<br>2/4   | MERO<br>1    | COL<br>2    | P/T4<br>32/4  | CZA<br>2/4   | MERO<br>1    | COL<br>2    | P/T4<br>32/4  | CZA<br>2/4   | MERO<br>1    |
| E | COL<br>4    | C/T<br>0.25/4 | CZA<br>4/4   | MERO<br>2    | COL<br>4    | C/T<br>0.25/4 | CZA<br>4/4   | MERO<br>2    | COL<br>4    | C/T<br>0.25/4 | CZA<br>4/4   | MERO<br>2    |
| F | COL<br>8    | C/T<br>0.5/4  | CZA<br>8/4   | MERO<br>4    | COL<br>8    | C/T<br>0.5/4  | CZA<br>8/4   | MERO<br>4    | COL<br>8    | C/T<br>0.5/4  | CZA<br>8/4   | MERO<br>4    |
| G | P/T4<br>1/4 | C/T<br>1/4    | CZA<br>16/4  | MERO<br>8    | P/T4<br>1/4 | C/T<br>1/4    | CZA<br>16/4  | MERO<br>8    | P/T4<br>1/4 | C/T<br>1/4    | CZA<br>16/4  | MERO<br>8    |
| H | P/T4<br>2/4 | C/T<br>2/4    | MERO<br>0.12 | MERO<br>16   | P/T4<br>2/4 | C/T<br>2/4    | MERO<br>0.12 | MERO<br>16   | P/T4<br>2/4 | C/T<br>2/4    | MERO<br>0.12 | MERO<br>16   |

## Antimicrobics

|             |                                      |
|-------------|--------------------------------------|
| <b>COL</b>  | Colistin                             |
| <b>P/T4</b> | Piperacillin / Tazobactam constant 4 |
| <b>C/T</b>  | Ceftolozane / Tazobactam 4           |
| <b>CZA</b>  | Ceftazidime / Avibactam              |
| <b>MERO</b> | Meropenem                            |
| <b>POS</b>  | Positive control                     |

\*For Proteus spp. \*\*For Enterobacteriaceae and Non-Enterobacteriaceae. \*\*\*For aid in detection of resistance mechanisms for Enterobacteriaceae and Non-Enterobacteriaceae.

# Sensititre™ Gram Negative IND1MDRF Plate

| Intended use                         | Organism                              | Use and methodology | Inoculum preparation                                                | Broth Type                               | Read method                                                                                      | Pack Size   |
|--------------------------------------|---------------------------------------|---------------------|---------------------------------------------------------------------|------------------------------------------|--------------------------------------------------------------------------------------------------|-------------|
| Antimicrobial Susceptibility Testing | Non-fastidious Gram Negative isolates | RUO (CLSI)          | Sensititre™ Sterile Water (T3339)<br>0.5 McFarland Standard (E1041) | Sensititre™ Mueller Hinton Broth (T3462) | Sensititre™ ARIS HIQ<br>Sensititre™ OptiRead<br>Sensititre™ Vision<br>Sensititre™ Manual Viewbox | 10 Isolates |

Put 3-5 colonies into H<sub>2</sub>O to measure a 0.5 McFarland using the Nephelometer, mix 1 µL\*, 10 µL\*\*, or 30 µL\*\*\* of suspension into MHb

Inoculate plate with 50 µL volume per well of the suspension using the Sensititre™ AIM or Multi-Channel Pipette

Seal Sensititre™ plate and incubate at 34-36°C in a non-CO<sub>2</sub> incubator or Sensititre™ ARIS HIQ for 18-24 hours

Read automatically with ARIS HIQ or OptiRead, semi-automatically with Vision or manually with Manual Viewbox

|   | 1             | 2             | 3             | 4            | 5          | 6             | 7            | 8           | 9           | 10          | 11         | 12          |
|---|---------------|---------------|---------------|--------------|------------|---------------|--------------|-------------|-------------|-------------|------------|-------------|
| A | MERO<br>0.05  | MERO<br>1     | MERO<br>2     | MERO<br>4    | MERO<br>8  | MERO<br>16    | PLZ<br>1     | PLZ<br>2    | PLZ<br>4    | PLZ<br>8    |            |             |
| B | MEV<br>0.06/8 | MEV<br>0.12/8 | MEV<br>0.25/8 | MEV<br>0.5/8 | MEV<br>1/8 | MEV<br>2/8    | MEV<br>4/8   | MEV<br>8/8  | MEV<br>16/8 |             |            |             |
| C | IMI<br>0.5    | IMI<br>1      | IMI<br>2      | IMI<br>4     | IMI<br>8   | IMR<br>0.25/4 | IMR<br>0.5/4 | IMR<br>1/4  | IMR<br>2/4  | IMR<br>4/4  | IMR<br>8/4 | IMR<br>16/4 |
| D | DLX<br>0.12   | DLX<br>0.25   | DLX<br>0.5    | DLX<br>1     | DLX<br>2   | DLX<br>4      | FOS+<br>16   | FOS+<br>32  | FOS+<br>64  | FOS+<br>128 |            |             |
| E | ERV<br>0.03   | ERV<br>0.06   | ERV<br>0.12   | ERV<br>0.25  | ERV<br>0.5 | ERV<br>1      | ERV<br>2     | ERV<br>4    | ERV<br>8    |             |            |             |
| F | C/T<br>0.25/4 | C/T<br>0.5/4  | C/T<br>1/4    | C/T<br>2/4   | C/T<br>4/4 | C/T<br>8/4    | C/T<br>16/4  | C/T<br>32/4 | COL<br>1    | COL<br>2    | COL<br>4   | COL<br>8    |
| G | CZA<br>0.12/4 | CZA<br>0.25/4 | CZA<br>0.5/4  | CZA<br>1/4   | CZA<br>2/4 | CZA<br>4/4    | CZA<br>8/4   | CZA<br>16/4 | CZA<br>32/4 |             |            | POS         |
| H | SUD<br>0.12/4 | SUD<br>0.25/4 | SUD<br>0.5/4  | SUD<br>1/4   | SUD<br>2/4 | SUD<br>4/4    | SUD<br>8/4   | SUD<br>16/4 |             |             | POS        | POS         |


## Antimicrobics

|             |                                     |
|-------------|-------------------------------------|
| <b>CZA</b>  | Ceftazidime / Avibactam Constant 4  |
| <b>C/T</b>  | Ceftolozane / Tazobactam Constant 4 |
| <b>COL</b>  | Colistin                            |
| <b>DLX</b>  | Delafloxacin                        |
| <b>ERV</b>  | Eravacycline                        |
| <b>FOS+</b> | Fosfomycin+glucose-6-phosphate      |
| <b>IMI</b>  | Imipenem                            |
| <b>IMR</b>  | Imipenem / Relebactam Constant 4    |
| <b>MERO</b> | Meropenem                           |
| <b>MEV</b>  | Meropenem / Vaborbactam Constant 8  |
| <b>PLZ</b>  | Plazomicin sulfate                  |
| <b>SUD</b>  | Sulbactam / Durlabactam Constant 4  |
| <b>POS</b>  | Positive Control                    |

\*For Proteus spp. \*\*For Enterobacteriaceae and Non-Enterobacteriaceae. \*\*\*For aid in detection of resistance mechanisms for Enterobacteriaceae and Non-Enterobacteriaceae.

# Sensititre™ Gram Negative GNX3F Plate

for second-line testing including Colistin



Put 3-5 colonies into H<sub>2</sub>O to measure a 0.5 McFarland using the Nephelometer, Mix 1 µl\*, 10 µl\*\*, or 30 µl\*\*\* of suspension into MHB

Inoculate plate with 50 µl volume per well of the suspension using the Sensititre™ AIM or Multi-Channel Pipette

Seal Sensititre™ plate and incubate at 34-36°C in a non-CO<sub>2</sub> incubator or Sensititre™ ARIS HIQ for 18-24 hours

Automatically read with ARIS HIQ or OptiRead or manually read with Sensititre™ Vision or Sensititre™ Manual Viewbox

|   | 1         | 2         | 3           | 4           | 5              | 6         | 7         | 8            | 9            | 10           | 11            | 12            |
|---|-----------|-----------|-------------|-------------|----------------|-----------|-----------|--------------|--------------|--------------|---------------|---------------|
| A | AMI<br>4  | GEN<br>1  | TOB<br>1    | CIP<br>0.06 | SXT<br>0.5/9.5 | AZT<br>2  | FEP<br>2  | COL<br>0.25  | COL<br>0.5   | COL<br>1     | COL<br>2      | COL<br>4      |
| B | AMI<br>8  | GEN<br>2  | TOB<br>2    | CIP<br>0.12 | SXT<br>1/19    | AZT<br>4  | FEP<br>4  | POL<br>0.25  | POL<br>0.5   | POL<br>1     | POL<br>2      | POL<br>4      |
| C | AMI<br>16 | GEN<br>4  | TOB<br>4    | CIP<br>0.25 | SXT<br>2/38    | AZT<br>8  | FEP<br>8  | TAZ<br>1     | TAZ<br>2     | TAZ<br>4     | TAZ<br>8      | TAZ<br>16     |
| D | AMI<br>32 | GEN<br>8  | TOB<br>8    | CIP<br>0.5  | SXT<br>4/76    | AZT<br>16 | FEP<br>16 | FOT<br>2     | FOT<br>4     | FOT<br>8     | FOT<br>16     | FOT<br>32     |
| E | DOX<br>2  | MIN<br>2  | TGC<br>0.25 | CIP<br>1    | LEVO<br>1      | IMI<br>1  | MERO<br>1 | A/S2<br>4/2  | A/S2<br>8/4  | A/S2<br>16/8 | A/S2<br>32/16 | A/S2<br>64/32 |
| F | DOX<br>4  | MIN<br>4  | TGC<br>0.5  | CIP<br>2    | LEVO<br>2      | IMI<br>2  | MERO<br>2 | DOR<br>0.5   | DOR<br>1     | DOR<br>2     | DOR<br>4      | POS           |
| G | DOX<br>8  | MIN<br>8  | TGC<br>1    | TGC<br>4    | LEVO<br>4      | IMI<br>4  | MERO<br>4 | P/T4<br>8/4  | P/T4<br>16/4 | P/T4<br>32/4 | P/T4<br>64/4  | POS           |
| H | DOX<br>16 | MIN<br>16 | TGC<br>2    | TGC<br>8    | LEVO<br>8      | IMI<br>8  | MERO<br>8 | TIM2<br>16/2 | TIM2<br>32/2 | TIM2<br>64/2 | TIM2<br>128/2 | POS           |

## Antimicrobics

|      |                                          |
|------|------------------------------------------|
| AMI  | Amikacin                                 |
| DOX  | Doxycycline                              |
| GEN  | Gentamicin                               |
| MIN  | Minocycline                              |
| TOB  | Tobramycin                               |
| TGC  | Tigecycline                              |
| CIP  | Ciprofloxacin                            |
| SXT  | Trimethoprim / Sulfamethoxazole          |
| LEVO | Levofloxacin                             |
| AZT  | Aztreonam                                |
| IMI  | Imipenem                                 |
| FEP  | Cefepime                                 |
| MERO | Meropenem                                |
| COL  | Colistin                                 |
| POL  | Polymyxin B                              |
| TAZ  | Ceftazidime                              |
| FOT  | Cefotaxime                               |
| A/S2 | Ampicillin / sulbactam 2:1 ratio         |
| DOR  | Doripenem                                |
| P/T4 | Piperacillin / Tazobactam constant 4     |
| TIM2 | Ticarcillin / clavulanic acid constant 2 |
| POS  | Positive control                         |

\*For Proteus spp. \*\*For Enterobacteriaceae and Non-Enterobacteriaceae. \*\*\*For aid in detection of resistance mechanisms for Enterobacteriaceae and Non-Enterobacteriaceae.

# Sensititre™ Gram Negative GNX2F Plate

| Intended use                         | Organism                              | Use and methodology | Inoculum preparation                                                | Broth Type                               | Read method                                                                                      | Pack Size   |
|--------------------------------------|---------------------------------------|---------------------|---------------------------------------------------------------------|------------------------------------------|--------------------------------------------------------------------------------------------------|-------------|
| Antimicrobial Susceptibility Testing | Non-fastidious Gram Negative isolates | RUO (CLSI)          | Sensititre™ Sterile Water (T3339)<br>0.5 McFarland Standard (E1041) | Sensititre™ Mueller Hinton Broth (T3462) | Sensititre™ ARIS HiQ<br>Sensititre™ OptiRead<br>Sensititre™ Vision<br>Sensititre™ Manual Viewbox | 10 Isolates |

Put 3-5 colonies into H<sub>2</sub>O to measure a 0.5 McFarland using the Nephelometer, Mix 1 µl\*, 10 µl\*\*, or 30 µl\*\*\* of suspension into MHb

Inoculate plate with 50 µl volume per well of the suspension using the Sensititre™ AIM or Multi-Channel Pipette

Seal Sensititre™ plate and incubate at 34-36°C in a non-CO<sub>2</sub> incubator or Sensititre™ ARIS HiQ for 18-24 hours

Automatically read with ARIS HiQ or OptiRead or manually read with Sensititre™ Vision or Sensititre™ Manual Viewbox

|   | 1             | 2            | 3              | 4         | 5         | 6           | 7           | 8           | 9           | 10         | 11        | 12        |
|---|---------------|--------------|----------------|-----------|-----------|-------------|-------------|-------------|-------------|------------|-----------|-----------|
| A | AMI<br>4      | AZT<br>2     | SXT<br>0.5/9.5 | FEP<br>2  | LEVO<br>1 | CIP<br>0.25 | MERO<br>1   | DOR<br>0.12 | DOR<br>0.25 | DOR<br>0.5 | DOR<br>1  | DOR<br>2  |
| B | AMI<br>8      | AZT<br>4     | SXT<br>1/19    | FEP<br>4  | LEVO<br>2 | CIP<br>0.5  | MERO<br>2   | COL<br>0.25 | COL<br>0.5  | COL<br>1   | COL<br>2  | COL<br>4  |
| C | AMI<br>16     | AZT<br>8     | SXT<br>2/38    | FEP<br>8  | LEVO<br>4 | CIP<br>1    | MERO<br>4   | POL<br>0.25 | POL<br>0.5  | POL<br>1   | POL<br>2  | POL<br>4  |
| D | AMI<br>32     | AZT<br>16    | SXT<br>4/76    | FEP<br>16 | LEVO<br>8 | CIP<br>2    | MERO<br>8   | TAZ<br>1    | TAZ<br>2    | TAZ<br>4   | TAZ<br>8  | TAZ<br>16 |
| E | TIM2<br>16/2  | P/T4<br>8/4  | GEN<br>1       | TOB<br>1  | DOX<br>2  | MIN<br>2    | FOT<br>1    | FOT<br>2    | FOT<br>4    | FOT<br>8   | FOT<br>16 | FOT<br>32 |
| F | TIM2<br>32/2  | P/T4<br>16/4 | GEN<br>2       | TOB<br>2  | DOX<br>4  | MIN<br>4    | TGC<br>0.25 | TGC<br>0.5  | TGC<br>1    | TGC<br>2   | TGC<br>4  | TGC<br>8  |
| G | TIM2<br>64/2  | P/T4<br>32/4 | GEN<br>4       | TOB<br>4  | DOX<br>8  | MIN<br>8    | ETP<br>0.25 | ETP<br>0.5  | ETP<br>1    | ETP<br>2   | ETP<br>4  | POS       |
| H | TIM2<br>128/2 | P/T4<br>64/4 | GEN<br>8       | TOB<br>8  | DOX<br>16 | MIN<br>16   | IMI<br>1    | IMI<br>2    | IMI<br>4    | IMI<br>8   | POS       | POS       |

## Antimicrobics

|             |                                          |
|-------------|------------------------------------------|
| <b>AMI</b>  | Amikacin                                 |
| <b>TIM2</b> | Ticarcillin / Clavulanic acid constant 2 |
| <b>AZT</b>  | Aztreonam                                |
| <b>P/T4</b> | Piperacillin / Tazobactam constant 4     |
| <b>SXT</b>  | Trimethoprim / Sulfamethoxazole          |
| <b>GEN</b>  | Gentamicin                               |
| <b>FEP</b>  | Cefepime                                 |
| <b>TOB</b>  | Tobramycin                               |
| <b>LEVO</b> | Levofloxacin                             |
| <b>DOX</b>  | Doxycycline                              |
| <b>CIP</b>  | Ciprofloxacin                            |
| <b>MIN</b>  | Minocycline                              |
| <b>MERO</b> | Meropenem                                |
| <b>FOT</b>  | Cefotaxime                               |
| <b>TGC</b>  | Tigecycline                              |
| <b>ETP</b>  | Ertapenem                                |
| <b>IMI</b>  | Imipenem                                 |
| <b>DOR</b>  | Doripenem                                |
| <b>COL</b>  | Colistin                                 |
| <b>POL</b>  | Polymyxin B                              |
| <b>TAZ</b>  | Ceftazidime                              |
| <b>POS</b>  | Positive control                         |

\*For Proteus spp. \*\*For Enterobacteriaceae and Non-Enterobacteriaceae. \*\*\*For aid in detection of resistance mechanisms for Enterobacteriaceae and Non-Enterobacteriaceae.

# Sensititre™ YeastOne YO10 Plate

with anidulafungin and micafungin

| Intended use                         | Organism                                                                                                                                                           | Use and methodology | Inoculum preparation                                                | Broth Type                         | Read method                                      | Pack Size   |
|--------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|---------------------------------------------------------------------|------------------------------------|--------------------------------------------------|-------------|
| Antimicrobial Susceptibility Testing | Non-fastidious yeast including <i>Candida</i> species, <i>Cryptococcus</i> species, <i>Aspergillus</i> species and miscellaneous other rapid growing yeast species | IVD (CLSI)          | Sensititre™ Sterile Water (T3339)<br>0.5 McFarland Standard (E1041) | Sensititre™ YeastOne Broth (Y3462) | Sensititre™ Vizion<br>Sensititre™ Manual Viewbox | 10 Isolates |

Put 3-5 colonies into H<sub>2</sub>O to measure a 0.5 McFarland using the Nephelometer, mix 20 µL of suspension into Sensititre™ YeastOne Broth

Inoculate plate with 100 µL volume per well of the suspension using the Sensititre™ AIM or Multi-Channel Pipette

Seal Sensititre™ plate and incubate at 35°C in a non-CO<sub>2</sub> air incubator for 24 to 25\* hours

Read semi-automatically with Vizion or manually with Manual Viewbox

|   | 1            | 2            | 3           | 4           | 5           | 6           | 7          | 8        | 9        | 10       | 11        | 12         |
|---|--------------|--------------|-------------|-------------|-------------|-------------|------------|----------|----------|----------|-----------|------------|
| A | POS          | AND<br>0.015 | AND<br>0.03 | AND<br>0.06 | AND<br>0.12 | AND<br>0.25 | AND<br>0.5 | AND<br>1 | AND<br>2 | AND<br>4 | AND<br>8  | AB<br>0.12 |
| B | MF<br>0.008  | MF<br>0.015  | MF<br>0.03  | MF<br>0.06  | MF<br>0.12  | MF<br>0.25  | MF<br>0.5  | MF<br>1  | MF<br>2  | MF<br>4  | MF<br>8   | AB<br>0.25 |
| C | CAS<br>0.008 | CAS<br>0.015 | CAS<br>0.03 | CAS<br>0.06 | CAS<br>0.12 | CAS<br>0.25 | CAS<br>0.5 | CAS<br>1 | CAS<br>2 | CAS<br>4 | CAS<br>8  | AB<br>0.5  |
| D | FC<br>0.06   | FC<br>0.12   | FC<br>0.25  | FC<br>0.5   | FC<br>1     | FC<br>2     | FC<br>4    | FC<br>8  | FC<br>16 | FC<br>32 | FC<br>64  | AB<br>1    |
| E | PZ<br>0.008  | PZ<br>0.015  | PZ<br>0.03  | PZ<br>0.06  | PZ<br>0.12  | PZ<br>0.25  | PZ<br>0.5  | PZ<br>1  | PZ<br>2  | PZ<br>4  | PZ<br>8   | AB<br>2    |
| F | VOR<br>0.008 | VOR<br>0.015 | VOR<br>0.03 | VOR<br>0.06 | VOR<br>0.12 | VOR<br>0.25 | VOR<br>0.5 | VOR<br>1 | VOR<br>2 | VOR<br>4 | VOR<br>8  | AB<br>4    |
| G | IZ<br>0.015  | IZ<br>0.03   | IZ<br>0.06  | IZ<br>0.12  | IZ<br>0.25  | IZ<br>0.5   | IZ<br>1    | IZ<br>2  | IZ<br>4  | IZ<br>8  | IZ<br>16  | AB<br>8    |
| H | FZ<br>0.12   | FZ<br>0.25   | FZ<br>0.5   | FZ<br>1     | FZ<br>2     | FZ<br>4     | FZ<br>8    | FZ<br>16 | FZ<br>32 | FZ<br>64 | FZ<br>128 | FZ<br>256  |

## Antimicrobics

|     |                  |
|-----|------------------|
| FC  | 5-Flucytosine    |
| AB  | Amphotericin B   |
| AND | Anidulafungin    |
| CAS | Caspofungin      |
| FZ  | Fluconazole      |
| IZ  | Itraconazole     |
| MF  | Micafungin       |
| PZ  | Posaconazole     |
| VOR | Voriconazole     |
| POS | Positive control |

\*24-25 hours incubation for *Candida* species only. Please refer the instruction manual for other species.

# Sensititre™ *Mycobacterium tuberculosis* MYCOTB Plate

| Intended use                         | Organism                                   | Use and methodology | Inoculum preparation                                                                | Broth Type                                    | Read method                                      | Pack Size   |
|--------------------------------------|--------------------------------------------|---------------------|-------------------------------------------------------------------------------------|-----------------------------------------------|--------------------------------------------------|-------------|
| Antimicrobial Susceptibility Testing | <i>Mycobacterium tuberculosis</i> isolates | RUO                 | Sensititre™ Saline Tween with Glass Beads (T3490)<br>0.5 McFarland Standard (E1041) | Sensititre™ Middlebrook 7H9 with OADC (T3440) | Sensititre™ Vision<br>Sensititre™ Manual Viewbox | 10 Isolates |

Put 3-5 colonies into Saline Tween with Glass Beads to reach 0.5 McFarland Standard, Mix 100  $\mu$ L into Sensititre™ Middlebrook 7H9 with OADC

Inoculate plate with 100  $\mu$ L volume per well of the suspension using the Sensititre™ AIM or Multi-Channel Pipette

Seal Sensititre™ plate and incubate at 35-37°C in a non-CO<sub>2</sub> incubator for 10-21 days


Read semi-automatically with Vision or manually with Manual Viewbox

|   | 1        | 2        | 3        | 4        | 5        | 6        | 7       | 8       | 9       | 10       | 11      | 12      |
|---|----------|----------|----------|----------|----------|----------|---------|---------|---------|----------|---------|---------|
| A | OFL 32   | MXF 8    | RIF 16   | AMI 16   | STR 32   | RFB 16   | PAS 64  | ETH 40  | CYC 256 | INH 4    | KAN 40  | EMB 32  |
| B | OFL 16   | MXF 4    | RIF 8    | AMI 8    | STR 16   | RFB 8    | PAS 32  | ETH 20  | CYC 128 | INH 2    | KAN 20  | EMB 16  |
| C | OFL 8    | MXF 2    | RIF 4    | AMI 4    | STR 8    | RFB 4    | PAS 16  | ETH 10  | CYC 64  | INH 1    | KAN 10  | EMB 8   |
| D | OFL 4    | MXF 1    | RIF 2    | AMI 2    | STR 4    | RFB 2    | PAS 8   | ETH 5   | CYC 32  | INH 0.5  | KAN 5   | EMB 4   |
| E | OFL 2    | MXF 0.5  | RIF 1    | AMI 1    | STR 2    | RFB 1    | PAS 4   | ETH 2.5 | CYC 16  | INH 0.25 | KAN 2.5 | EMB 2   |
| F | OFL 1    | MXF 0.25 | RIF 0.5  | AMI 0.5  | STR 1    | RFB 0.5  | PAS 2   | ETH 1.2 | CYC 8   | INH 0.12 | KAN 1.2 | EMB 1   |
| G | OFL 0.5  | MXF 0.12 | RIF 0.25 | AMI 0.25 | STR 0.5  | RFB 0.25 | PAS 1   | ETH 0.6 | CYC 4   | INH 0.06 | KAN 0.6 | EMB 0.5 |
| H | OFL 0.25 | MXF 0.06 | RIF 0.12 | AMI 0.12 | STR 0.25 | RFB 0.12 | PAS 0.5 | ETH 0.3 | CYC 2   | INH 0.03 | POS     | POS     |

## Antimicrobics

|            |                          |
|------------|--------------------------|
| <b>AMI</b> | Amikacin                 |
| <b>CYC</b> | Cycloserine              |
| <b>EMB</b> | Ethambutol               |
| <b>ETH</b> | Ethionamide              |
| <b>INH</b> | Isoniazid                |
| <b>KAN</b> | Kanamycin                |
| <b>MXF</b> | Moxifloxacin             |
| <b>OFL</b> | Ofloxacin                |
| <b>PAS</b> | Para-aminosalicylic acid |
| <b>RFB</b> | Rifabutin                |
| <b>RIF</b> | Rifampin                 |
| <b>STR</b> | Streptomycin             |
| <b>POS</b> | Positive control         |

# Sensititre™ *Mycobacterium tuberculosis* IND1MBT Plate



Put 3-5 colonies into Saline Tween with Glass Beads to reach 0.5 McFarland Standard, Mix 100 µL into Sensititre™ Middlebrook 7H9 with OADC

Inoculate plate with 100 µL volume per well of the suspension using the Sensititre™ AIM or Multi-Channel Pipette

Seal Sensititre™ plate and incubate at 35-37°C in a non-CO<sub>2</sub> incubator for 10-21 days

Read semi-automatically with Vision or manually with Manual Viewbox

|   | 1           | 2           | 3           | 4          | 5           | 6           | 7           | 8           | 9           | 10           | 11           | 12       |
|---|-------------|-------------|-------------|------------|-------------|-------------|-------------|-------------|-------------|--------------|--------------|----------|
| A | RIF<br>0.06 | INH<br>0.03 | MXF<br>0.06 | ETH<br>0.3 | LZD<br>0.03 | CFZ<br>0.06 | PAS<br>0.12 | CYC<br>0.06 | EMB<br>0.25 | LEVO<br>0.12 | BDQ<br>0.015 | BDQ<br>2 |
| B | RIF<br>0.12 | INH<br>0.06 | MXF<br>0.12 | ETH<br>0.6 | LZD<br>0.06 | CFZ<br>0.12 | PAS<br>0.25 | CYC<br>0.12 | EMB<br>0.5  | LEVO<br>0.25 | BDQ<br>0.03  | BDQ<br>4 |
| C | RIF<br>0.25 | INH<br>0.12 | MXF<br>0.25 | ETH<br>1.2 | LZD<br>0.12 | CFZ<br>0.25 | PAS<br>0.5  | CYC<br>0.25 | EMB<br>1    | LEVO<br>0.5  | BDQ<br>0.06  |          |
| D | RIF<br>0.5  | INH<br>0.25 | MXF<br>0.5  | ETH<br>2.5 | LZD<br>0.25 | CFZ<br>0.5  | PAS<br>1    | CYC<br>0.5  | EMB<br>2    | LEVO<br>1    | BDQ<br>0.12  |          |
| E | RIF<br>1    | INH<br>0.5  | MXF<br>1    | ETH<br>5   | LZD<br>5    | CFZ<br>1    | PAS<br>2    | CYC<br>1    | EMB<br>4    | LEVO<br>2    | BDQ<br>0.25  |          |
| F | RIF<br>2    | INH<br>1    | MXF<br>2    | ETH<br>10  | LZD<br>1    | CFZ<br>2    | PAS<br>4    | CYC<br>2    | EMB<br>8    | LEVO<br>4    | BDQ<br>0.5   |          |
| G | RIF<br>4    | INH<br>2    | MXF<br>4    | ETH<br>20  | LZD<br>2    | CFZ<br>4    | PAS<br>8    | CYC<br>4    | EMB<br>16   | LEVO<br>8    | BDQ<br>1     |          |
| H | RIF<br>8    | INH<br>4    | MXF<br>8    | AMI<br>40  | AMI<br>0.12 | AMI<br>0.25 | AMI<br>0.5  | AMI<br>1    | AMI<br>2    | AMI<br>4     | AMI<br>8     | POS      |

## Antimicrobics

|             |                          |
|-------------|--------------------------|
| <b>RIF</b>  | Rifampin                 |
| <b>INH</b>  | Isoniazid                |
| <b>MXF</b>  | Moxifloxacin             |
| <b>ETH</b>  | Ethionamide              |
| <b>LZD</b>  | Linezolid                |
| <b>AMI</b>  | Amikacin                 |
| <b>CFZ</b>  | Clofazimine              |
| <b>PAS</b>  | Para-aminosalicylic acid |
| <b>CYC</b>  | Cycloserine              |
| <b>EMB</b>  | Ethambutol               |
| <b>LEVO</b> | Levofloxacin             |
| <b>BDQ</b>  | Bedaquiline              |
| <b>POS</b>  | Positive Control         |

# Sensititre™ Rapid Growing Mycobacteria RAPMYCOI Plate

| Intended use                         | Organism                                                                                                   | Use and methodology | Inoculum preparation                                                 | Broth Type                       | Read method                                      | Pack Size   |
|--------------------------------------|------------------------------------------------------------------------------------------------------------|---------------------|----------------------------------------------------------------------|----------------------------------|--------------------------------------------------|-------------|
| Antimicrobial Susceptibility Testing | Rapidly growing <i>Mycobacteria</i> , <i>Nocardia</i> spp. and other aerobic <i>Actinomycetes</i> Isolates | IVD (CLSI)          | Sensititre™ Sterile Water (T3339)*<br>0.5 McFarland Standard (E1041) | Sensititre™ MHB with TES (T3462) | Sensititre™ Vision<br>Sensititre™ Manual Viewbox | 10 Isolates |

Put 3-5 colonies into H<sub>2</sub>O to measure a 0.5 McFarland using the Nephelometer, mix 50 µL of suspension into Sensititre™ MHB with TES

Inoculate plate with 100 µL volume per well of the suspension using the Sensititre™ AIM or Multi-Channel Pipette

Seal Sensititre™ plate and incubate at 30°C in a non-CO<sub>2</sub> air incubator for 72 hours\*\*

Read semi-automatically with Vision or manually with Manual Viewbox

|   | 1                | 2              | 3           | 4           | 5           | 6            | 7           | 8           | 9           | 10           | 11            | 12            |
|---|------------------|----------------|-------------|-------------|-------------|--------------|-------------|-------------|-------------|--------------|---------------|---------------|
| A | SXT<br>0.25/4.75 | SXT<br>0.5/9.5 | SXT<br>1/19 | SXT<br>2/38 | SXT<br>4/76 | SXT<br>8/152 | LZD<br>1    | LZD<br>2    | LZD<br>4    | LZD<br>8     | LZD<br>16     | LZD<br>32     |
| B | CIP<br>0.12      | CIP<br>0.25    | CIP<br>0.5  | CIP<br>1    | CIP<br>2    | CIP<br>4     | IMI<br>2    | IMI<br>4    | IMI<br>8    | IMI<br>16    | IMI<br>32     | IMI<br>64     |
| C | MXF<br>0.25      | MXF<br>0.5     | MXF<br>1    | MXF<br>2    | MXF<br>4    | MXF<br>8     | FEP<br>1    | FEP<br>2    | FEP<br>4    | FEP<br>8     | FEP<br>16     | FEP<br>32     |
| D | FOX<br>4         | FOX<br>8       | FOX<br>16   | FOX<br>32   | FOX<br>64   | FOX<br>128   | AUG2<br>2/1 | AUG2<br>4/2 | AUG2<br>8/4 | AUG2<br>16/8 | AUG2<br>32/16 | AUG2<br>64/32 |
| E | AMI<br>1         | AMI<br>2       | AMI<br>4    | AMI<br>8    | AMI<br>16   | AMI<br>32    | AMI<br>64   | AXO<br>4    | AXO<br>8    | AXO<br>16    | AXO<br>32     | AXO<br>64     |
| F | DOX<br>0.12      | DOX<br>0.25    | DOX<br>0.5  | DOX<br>1    | DOX<br>2    | DOX<br>4     | DOX<br>8    | DOX<br>16   | MIN<br>1    | MIN<br>2     | MIN<br>4      | MIN<br>8      |
| G | TGC<br>0.015     | TGC<br>0.03    | TGC<br>0.06 | TGC<br>0.12 | TGC<br>0.25 | TGC<br>0.5   | TGC<br>1    | TGC<br>2    | TGC<br>4    | TOB<br>1     | TOB<br>2      | TOB<br>4      |
| H | CLA<br>0.06      | CLA<br>0.12    | CLA<br>0.25 | CLA<br>0.5  | CLA<br>1    | CLA<br>2     | CLA<br>4    | CLA<br>8    | CLA<br>16   | TOB<br>8     | TOB<br>16     | POS           |

## Antimicrobics

|      |                                       |
|------|---------------------------------------|
| AMI  | Amikacin                              |
| AUG2 | Amoxicillin/Clavulanic acid 2:1 ratio |
| FEP  | Cefepime                              |
| FOX  | Cefoxitin                             |
| AXO  | Ceftriaxone                           |
| CIP  | Ciprofloxacin                         |
| CLA  | Clarithromycin                        |
| DOX  | Doxycycline                           |
| IMI  | Imipenem                              |
| LZD  | Linezolid                             |
| MIN  | Minocycline                           |
| MXF  | Moxifloxacin                          |
| TGC  | Tigecycline                           |
| TOB  | Tobramycin                            |
| SXT  | Trimethoprim/Sulfamethoxazole         |
| POS  | Positive control                      |

\*For Hard colonies, vortex with demineralized water with glass beads (T3492) to make a homogenous suspension.

\*\*Refer instruction manual for incubation specification for different isolates.

# Sensititre™ Rapid Growing Mycobacteria RAPMYCO2 Plate

| Intended use                         | Organism                                                                                                   | Use and methodology | Inoculum preparation                                                 | Broth Type                       | Read method                                      | Pack Size   |
|--------------------------------------|------------------------------------------------------------------------------------------------------------|---------------------|----------------------------------------------------------------------|----------------------------------|--------------------------------------------------|-------------|
| Antimicrobial Susceptibility Testing | Rapidly growing <i>Mycobacteria</i> , <i>Nocardia</i> spp. and other aerobic <i>Actinomycetes</i> Isolates | RUO                 | Sensititre™ Sterile Water (T3339)*<br>0.5 McFarland Standard (E1041) | Sensititre™ MHB with TES (T3462) | Sensititre™ Vision<br>Sensititre™ Manual Viewbox | 10 Isolates |

Put 3-5 colonies into H<sub>2</sub>O to measure a 0.5 McFarland using the Nephelometer, mix 50 µL of suspension into Sensititre™ MHB with TES

Inoculate plate with 100 µL volume per well of the suspension using the Sensititre™ AIM or Multi-Channel Pipette

Seal Sensititre™ plate and incubate at 30°C in a non-CO<sub>2</sub> air incubator for 72 hours\*\*

Read semi-automatically with Vision or manually with Manual Viewbox

|   | 1            | 2            | 3           | 4           | 5           | 6           | 7                | 8              | 9           | 10          | 11          | 12          |
|---|--------------|--------------|-------------|-------------|-------------|-------------|------------------|----------------|-------------|-------------|-------------|-------------|
| A | IMI<br>0.008 | IMI<br>0.015 | IMI<br>0.03 | IMI<br>0.06 | IMI<br>0.12 | IMI<br>0.25 | IMI<br>0.5       | IMI<br>1       | IMI<br>2    | IMI<br>4    | IMI<br>8    | IMI<br>16   |
| B | IMI<br>32    | TGC<br>0.03  | TGC<br>0.06 | TGC<br>0.12 | TGC<br>0.25 | TGC<br>0.5  | TGC<br>1         | TGC<br>2       | DOX<br>0.12 | DOX<br>0.25 | DOX<br>0.5  | DOX<br>1    |
| C | DOX<br>2     | DOX<br>4     | DOX<br>8    | AMI<br>1    | AMI<br>2    | AMI<br>4    | AMI<br>8         | AMI<br>16      | AMI<br>32   | AMI<br>64   | AMI<br>128  | AMI<br>256  |
| D | CIP<br>0.12  | CIP<br>0.25  | CIP<br>0.5  | CIP<br>1    | CIP<br>2    | CIP<br>4    | LZD<br>1         | LZD<br>2       | LZD<br>4    | LZD<br>8    | LZD<br>16   | LZD<br>32   |
| E | TOB<br>0.12  | TOB<br>0.25  | TOB<br>0.5  | TOB<br>1    | TOB<br>2    | TOB<br>4    | TOB<br>8         | TOB<br>16      | FOX<br>1    | FOX<br>2    | FOX<br>4    | FOX<br>8    |
| F | FOX<br>16    | FOX<br>32    | FOX<br>64   | FOX<br>128  | CFZ<br>0.03 | CFZ<br>0.06 | CFZ<br>0.12      | CFZ<br>0.25    | CFZ<br>0.5  | CFZ<br>1    | CFZ<br>2    | CFZ<br>4    |
| G | MXF<br>0.015 | MXF<br>0.03  | MXF<br>0.06 | MXF<br>0.12 | MXF<br>0.25 | MXF<br>0.5  | MXF<br>1         | MXF<br>2       | MXF<br>4    | CLA<br>0.06 | CLA<br>0.12 | CLA<br>0.25 |
| H | CLA<br>0.5   | CLA<br>1     | CLA<br>2    | CLA<br>4    | CLA<br>8    | CLA<br>16   | SXT<br>0.25/4/75 | SXT<br>0.5/9.5 | SXT<br>1/19 | SXT<br>2/38 | SXT<br>4/76 | POS         |

## Antimicrobics

|            |                                 |
|------------|---------------------------------|
| <b>IMI</b> | Imipenem                        |
| <b>DOX</b> | Doxycycline                     |
| <b>CIP</b> | Ciprofloxacin                   |
| <b>TOB</b> | Tobramycin                      |
| <b>FOX</b> | Cefoxitin                       |
| <b>MXF</b> | Moxifloxacin                    |
| <b>CLA</b> | Clarithromycin                  |
| <b>TGC</b> | Tigecycline                     |
| <b>AMI</b> | Amikacin                        |
| <b>CFZ</b> | Clofazimine                     |
| <b>LZD</b> | Linezolid                       |
| <b>SXT</b> | Trimethoprim / sulfamethoxazole |
| <b>POS</b> | Positive Control                |

\*For Hard colonies, vortex with demineralized water with glass beads (T3492) to make a homogenous suspension.

\*\*Refer instruction manual for incubation specification for different isolates.

# Sensititre™ Slow Growing Mycobacteria SLOMYCOI Plate

| Intended use                         | Organism                                                                                                                                                             | Use and methodology | Inoculum preparation                                                | Broth Type                        | Read method                                      | Pack Size   |
|--------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|---------------------------------------------------------------------|-----------------------------------|--------------------------------------------------|-------------|
| Antimicrobial Susceptibility Testing | Slowly growing non-tuberculosis mycobacteria (NTM), i.e. <i>Mycobacterium avium</i> complex, <i>Mycobacterium kansasii</i> and <i>Mycobacterium marinum</i> isolates | IVD (CLSI)          | Sensititre™ Sterile Water (T3339)<br>0.5 McFarland Standard (E1041) | Sensititre™ MHB with OADC (T8005) | Sensititre™ Vizion<br>Sensititre™ Manual Viewbox | 10 Isolates |

Put 3-5 colonies into H<sub>2</sub>O to measure a 0.5 McFarland using the Nephelometer, mix 50 µL of suspension into Sensititre™ MHB with OADC

Inoculate plate with 100 µL volume per well of the suspension using the Sensititre™ AIM or Multi-Channel Pipette

Seal Sensititre™ plate and incubate at 35°C in a non-CO<sub>2</sub> air incubator for 7-14 days

Read semi-automatically with Vizion or manually with Manual Viewbox

|   | 1           | 2           | 3           | 4           | 5           | 6                | 7         | 8         | 9           | 10         | 11          | 12         |
|---|-------------|-------------|-------------|-------------|-------------|------------------|-----------|-----------|-------------|------------|-------------|------------|
| A | CLA<br>0.06 | CLA<br>0.12 | CLA<br>0.25 | CLA<br>0.5  | CLA<br>1    | CLA<br>2         | CLA<br>4  | CLA<br>8  | CIP<br>16   | STR<br>64  | DOX<br>16   | ETH<br>20  |
| B | CLA<br>16   | CLA<br>32   | CLA<br>64   | MXF<br>8    | RIF<br>8    | SXT<br>8/152     | AMI<br>64 | LZD<br>64 | CIP<br>8    | STR<br>32  | DOX<br>8    | ETH<br>10  |
| C | RFB<br>8    | EMB<br>16   | INH<br>8    | MXF<br>4    | RIF<br>4    | SXT<br>4/76      | AMI<br>32 | LZD<br>32 | CIP<br>4    | STR<br>16  | DOX<br>4    | ETH<br>5   |
| D | RFB<br>4    | EMB<br>8    | INH<br>4    | MXF<br>2    | RIF<br>2    | SXT<br>2/38      | AMI<br>16 | LZD<br>16 | CIP<br>2    | STR<br>8   | DOX<br>2    | ETH<br>2.5 |
| E | RFB<br>2    | EMB<br>4    | INH<br>2    | MXF<br>1    | RIF<br>1    | SXT<br>1/19      | AMI<br>8  | LZD<br>8  | CIP<br>1    | STR<br>4   | DOX<br>1    | ETH<br>1.2 |
| F | RFB<br>1    | EMB<br>2    | INH<br>1    | MXF<br>0.5  | RIF<br>0.5  | SXT<br>0.5/9.5   | AMI<br>4  | LZD<br>4  | CIP<br>0.5  | STR<br>2   | DOX<br>0.5  | ETH<br>0.6 |
| G | RFB<br>0.5  | EMB<br>1    | INH<br>0.5  | MXF<br>0.25 | RIF<br>0.25 | SXT<br>0.25/4.75 | AMI<br>2  | LZD<br>2  | CIP<br>0.25 | STR<br>1   | DOX<br>0.25 | ETH<br>0.3 |
| H | RFB<br>0.25 | EMB<br>0.5  | INH<br>0.25 | MXF<br>0.12 | RIF<br>0.12 | SXT<br>0.12/2.38 | AMI<br>1  | LZD<br>1  | CIP<br>0.12 | STR<br>0.5 | DOX<br>0.12 | POS        |

## Antimicrobics

|            |                               |
|------------|-------------------------------|
| <b>AMI</b> | Amikacin                      |
| <b>CIP</b> | Ciprofloxacin                 |
| <b>CLA</b> | Clarithromycin                |
| <b>DOX</b> | Doxycycline                   |
| <b>EMB</b> | Ethambutol                    |
| <b>ETH</b> | Ethionamide                   |
| <b>INH</b> | Isoniazid                     |
| <b>LZD</b> | Linezolid                     |
| <b>MXF</b> | Moxifloxacin                  |
| <b>RFB</b> | Rifabutin                     |
| <b>RIF</b> | Rifampin                      |
| <b>STR</b> | Streptomycin                  |
| <b>SXT</b> | Trimethoprim/Sulfamethoxazole |
| <b>POS</b> | Positive control              |

# Sensititre™ Slow Growing Mycobacteria SLOMYCO2 Plate

| Intended use                         | Organism                                                                                                                                                             | Use and methodology | Inoculum preparation                                               | Broth Type                        | Read method                                      | Pack Size   |
|--------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|--------------------------------------------------------------------|-----------------------------------|--------------------------------------------------|-------------|
| Antimicrobial Susceptibility Testing | Slowly growing non-tuberculosis mycobacteria (NTM), i.e. <i>Mycobacterium avium</i> complex, <i>Mycobacterium kansasii</i> and <i>Mycobacterium marinum</i> isolates | IVD (CLSI)          | Sensititre Sterile Water (T3339)<br>0.5 McFarland Standard (E1041) | Sensititre™ MHB with OADC (T8005) | Sensititre™ Vision<br>Sensititre™ Manual Viewbox | 10 Isolates |

Put 3-5 colonies into H<sub>2</sub>O to measure a 0.5 McFarland using the Nephelometer, mix 50 µL of suspension into Sensititre™ MHB with OADC

Inoculate plate with 100 µL volume per well of the suspension using the Sensititre AIM or Multi-Channel Pipette

Seal Sensititre™ plate and incubate at 35°C in a non-CO<sub>2</sub> air incubator for 7-14 days

Read semi-automatically with Vision or manually with Manual Viewbox

|   | 1                | 2              | 3            | 4           | 5           | 6           | 7            | 8           | 9           | 10           | 11          | 12          |
|---|------------------|----------------|--------------|-------------|-------------|-------------|--------------|-------------|-------------|--------------|-------------|-------------|
| A | SXT<br>0.25/4.75 | SXT<br>0.5/9.5 | SXT<br>1/19  | SXT<br>2/38 | SXT<br>4/76 | DOX<br>0.12 | DOX<br>0.25  | DOX<br>0.5  | DOX<br>1    | DOX<br>2     | DOX<br>4    | DOX<br>8    |
| B | LZD<br>1         | LZD<br>2       | LZD<br>4     | LZD<br>8    | LZD<br>16   | LZD<br>32   | RFB<br>0.12  | RFB<br>0.25 | RFB<br>0.5  | RFB<br>1     | RFB<br>2    | RFB<br>4    |
| C | AMI<br>1         | AMI<br>2       | AMI<br>4     | AMI<br>8    | AMI<br>16   | AMI<br>32   | AMI<br>64    | AMI<br>128  | AMI<br>256  | MXF<br>0.015 | MXF<br>0.03 | MXF<br>0.06 |
| D | MXF<br>0.12      | MXF<br>0.25    | MXF<br>0.5   | MXF<br>1    | MXF<br>2    | MXF<br>4    | CFZ<br>0.015 | CFZ<br>0.03 | CFZ<br>0.06 | CFZ<br>0.12  | CFZ<br>0.25 | CFZ<br>0.5  |
| E | CFZ<br>1         | CFZ<br>2       | CFZ<br>4     | CIP<br>0.12 | CIP<br>0.25 | CIP<br>0.5  | CIP<br>1     | CIP<br>2    | CIP<br>4    | CIP<br>8     | STR<br>0.5  | STR<br>1    |
| F | STR<br>2         | STR<br>4       | STR<br>8     | STR<br>16   | STR<br>32   | CLA<br>0.06 | CLA<br>0.12  | CLA<br>0.25 | CLA<br>0.5  | CLA<br>1     | CLA<br>2    | CLA<br>4    |
| G | CLA<br>8         | CLA<br>16      | CLA<br>32    | CLA<br>64   | MIN<br>0.06 | MIN<br>0.12 | MIN<br>0.25  | MIN<br>0.5  | MIN<br>1    | MIN<br>2     | MIN<br>4    | MIN<br>8    |
| H | RIF<br>0.004     | RIF<br>0.008   | RIF<br>0.015 | RIF<br>0.03 | RIF<br>0.06 | RIF<br>0.12 | RIF<br>0.25  | RIF<br>0.5  | RIF<br>1    | RIF<br>2     | RIF<br>4    | POS         |

## Antimicrobics

|            |                                 |
|------------|---------------------------------|
| <b>SXT</b> | Trimethoprim / sulfamethoxazole |
| <b>LZD</b> | Linezolid                       |
| <b>AMI</b> | Amikacin                        |
| <b>MXF</b> | Moxifloxacin                    |
| <b>CFZ</b> | Clofazimine                     |
| <b>STR</b> | Streptomycin                    |
| <b>CLA</b> | Clarithromycin                  |
| <b>RIF</b> | Rifampin                        |
| <b>CIP</b> | Ciprofloxacin                   |
| <b>MIN</b> | Minocycline                     |
| <b>DOX</b> | Doxycycline                     |
| <b>RFB</b> | Rifabutin                       |
| <b>POS</b> | Positive Control                |

# Sensititre™ NOCARDIA Plate

| Intended use                         | Organism                                   | Use and methodology | Inoculum preparation                                                | Broth Type                       | Read method                                      | Pack Size   |
|--------------------------------------|--------------------------------------------|---------------------|---------------------------------------------------------------------|----------------------------------|--------------------------------------------------|-------------|
| Antimicrobial Susceptibility Testing | Susceptibility testing of Nocardia species | RUO                 | Sensititre™ Sterile Water (T3492)<br>0.5 McFarland Standard (E1041) | Sensititre™ MHB with TES (T3462) | Sensititre™ Vision<br>Sensititre™ Manual Viewbox | 10 Isolates |

Put 3-5 colonies into H<sub>2</sub>O to measure a 0.5 McFarland using the Nephelometer, mix 50 µL of suspension into Sensititre™ MHB with TES

Inoculate plate with 100 µL volume per well of the suspension using the Sensititre™ AIM or Multi-Channel Pipette

Seal Sensititre™ plate and incubate at 35°C in a non-CO<sub>2</sub> air incubator for 2-3 days

Read semi-automatically with Vision or manually with Manual Viewbox

|   | 1                | 2                | 3                | 4              | 5             | 6             | 7              | 8            | 9           | 10          | 11          | 12          |
|---|------------------|------------------|------------------|----------------|---------------|---------------|----------------|--------------|-------------|-------------|-------------|-------------|
| A | TOB<br>2         | TOB<br>4         | TOB<br>8         | TOB<br>16      | TOB<br>32     | TOB<br>64     | LZD<br>0.5     | LZD<br>1     | LZD<br>2    | LZD<br>4    | LZD<br>8    | LZD<br>16   |
| B | SXT<br>0.06/1.19 | SXT<br>0.12/2.38 | SXT<br>0.25/4.75 | SXT<br>0.5/9.5 | SXT<br>1/19   | SXT<br>2/38   | SXT<br>4/76    | SXT<br>8/152 | DOX<br>0.12 | DOX<br>0.25 | DOX<br>0.5  | DOX<br>1    |
| C | DOX<br>2         | DOX<br>4         | DOX<br>8         | DOX<br>16      | CIP<br>0.12   | CIP<br>0.25   | CIP<br>0.5     | CIP<br>1     | CIP<br>2    | CIP<br>4    | CIP<br>8    | CIP<br>16   |
| D | AUG2<br>2/1      | AUG2<br>4/2      | AUG2<br>8/4      | AUG2<br>16/8   | AUG2<br>32/16 | AUG2<br>64/32 | AUG2<br>128/64 | MXF<br>0.015 | MXF<br>0.03 | MXF<br>0.06 | MXF<br>0.12 | MXF<br>0.25 |
| E | MXF<br>0.5       | MXF<br>1         | MXF<br>2         | MXF<br>4       | MXF<br>8      | AMI<br>0.5    | AMI<br>1       | AMI<br>2     | AMI<br>4    | AMI<br>8    | AMI<br>16   | AMI<br>32   |
| F | IMI<br>0.12      | IMI<br>0.25      | IMI<br>0.5       | IMI<br>1       | IMI<br>2      | IMI<br>4      | IMI<br>8       | IMI<br>16    | IMI<br>32   | AXO<br>1    | AXO<br>2    | AXO<br>4    |
| G | AXO<br>8         | AXO<br>16        | AXO<br>32        | AXO<br>64      | CLA<br>0.03   | CLA<br>0.06   | CLA<br>0.12    | CLA<br>0.25  | CLA<br>0.5  | CLA<br>1    | CLA<br>2    | CLA<br>4    |
| H | CLA<br>8         | CLA<br>16        | MIN<br>0.06      | MIN<br>0.12    | MIN<br>0.25   | MIN<br>0.5    | MIN<br>1       | MIN<br>2     | MIN<br>4    | MIN<br>8    | MIN<br>16   | POS         |

## Antimicrobics

|             |                                         |
|-------------|-----------------------------------------|
| <b>TOM</b>  | Tobramycin                              |
| <b>SXT</b>  | Trimethoprim / sulfamethoxazole         |
| <b>DOX</b>  | Doxycycline                             |
| <b>AUG2</b> | Amoxicillin / clavulanic acid 2:1 ratio |
| <b>MXF</b>  | Moxifloxacin                            |
| <b>IMI</b>  | Imipenem                                |
| <b>AXO</b>  | Ceftriaxone                             |
| <b>CLA</b>  | Clarithromycin                          |
| <b>MIN</b>  | Minocycline                             |
| <b>CIP</b>  | Ciprofloxacin                           |
| <b>AMI</b>  | Amikacin                                |
| <b>LZD</b>  | Linezolid                               |
| <b>POS</b>  | Positive Control                        |

# Sensititre™ Anaerobe ANO2B Plate

| Intended use                         | Organism                                                | Use and methodology | Inoculum preparation                                                                                   | Broth Type                                                    | Read method                                      | Pack Size   |
|--------------------------------------|---------------------------------------------------------|---------------------|--------------------------------------------------------------------------------------------------------|---------------------------------------------------------------|--------------------------------------------------|-------------|
| Antimicrobial Susceptibility Testing | Anerobic organisms<br><i>Bacteroides fragilis</i> group | RUO                 | Sensititre™ Cation Adjusted Mueller Hinton Broth with TES (T3462-05)<br>0.5 McFarland Standard (E1041) | Sensititre™ Supplemented Brucella Broth for Anaerobes (T3451) | Sensititre™ Vision<br>Sensititre™ Manual Viewbox | 10 Isolates |

Put 3-5 colonies into 5ml MHB to measure a 0.5 McFarland using the Nephelometer, mix 100 µL of suspension into supplemented Brucella Broth

Inoculate plate with 100 µL volume per well of the suspension using the Sensititre™ AIM or Multi-Channel Pipette

Seal Sensititre™ plate and incubate at 34-36°C in an anaerobic atmosphere for 46-48 hours

Read semi-automatically with Vision or manually with Manual Viewbox

|   | 1                | 2             | 3           | 4           | 5           | 6            | 7                | 8             | 9            | 10            | 11          | 12           |
|---|------------------|---------------|-------------|-------------|-------------|--------------|------------------|---------------|--------------|---------------|-------------|--------------|
| A | A/S2<br>0.5/0.25 | A/S2<br>1/0.5 | A/S2<br>2/1 | A/S2<br>4/2 | A/S2<br>8/4 | A/S2<br>16/8 | AUG2<br>0.5/0.25 | AUG2<br>1/0.5 | AUG2<br>2/1  | AUG2<br>4/2   | AUG2<br>8/4 | AUG2<br>16/8 |
| B | TANS<br>4        | TANS<br>8     | TANS<br>16  | TANS<br>32  | TANS<br>64  | PEN<br>0.06  | PEN<br>0.12      | PEN<br>0.25   | PEN<br>0.5   | PEN<br>1      | PEN<br>2    | PEN<br>4     |
| C | IMI<br>0.12      | IMI<br>0.25   | IMI<br>0.5  | IMI<br>1    | IMI<br>2    | IMI<br>4     | IMI<br>8         | MERO<br>0.5   | MERO<br>1    | MERO<br>2     | MERO<br>4   | MERO<br>8    |
| D | CLI<br>0.25      | CLI<br>0.5    | CLI<br>1    | CLI<br>2    | CLI<br>4    | CLI<br>8     | FOX<br>1         | FOX<br>2      | FOX<br>4     | FOX<br>8      | FOX<br>16   | FOX<br>32    |
| E | MRD<br>0.5       | MRD<br>1      | MRD<br>2    | MRD<br>4    | MRD<br>8    | MRD<br>16    | CHL<br>2         | CHL<br>4      | CHL<br>8     | CHL<br>16     | CHL<br>32   | CHL<br>64    |
| F | AMP<br>0.5       | AMP<br>1      | AMP<br>2    | AMP<br>4    | AMP<br>8    | AMP<br>16    | PIP<br>4         | PIP<br>8      | PIP<br>16    | PIP<br>32     | PIP<br>64   | PIP<br>128   |
| G | TET<br>0.25      | TET<br>0.5    | TET<br>1    | TET<br>2    | TET<br>4    | TET<br>8     | MEZ<br>4         | MEZ<br>8      | MEZ<br>16    | MEZ<br>32     | MEZ<br>64   | MEZ<br>128   |
| H | P/T4<br>0.25/4   | P/T4<br>0.5/4 | P/T4<br>1/4 | P/T4<br>2/4 | P/T4<br>4/4 | P/T4<br>8/4  | P/T4<br>16/4     | P/T4<br>32/4  | P/T4<br>64/4 | P/T4<br>128/4 | POS         | POS          |

## Antimicrobics

|             |                                        |
|-------------|----------------------------------------|
| <b>AUG2</b> | Amoxicillin /Clavulanic acid 2:1 ratio |
| <b>AMP</b>  | Ampicillin                             |
| <b>A/S2</b> | Ampicillin/Sulbactam 2:1 ratio         |
| <b>TANS</b> | Cefotetan                              |
| <b>FOX</b>  | Cefoxitin                              |
| <b>CHL</b>  | Chloramphenicol                        |
| <b>CLI</b>  | Clindamycin                            |
| <b>IMI</b>  | Imipenem                               |
| <b>MERO</b> | Meropenem                              |
| <b>MRD</b>  | Metronidazole                          |
| <b>MEZ</b>  | Mezlocillin                            |
| <b>PEN</b>  | Penicillin                             |
| <b>PIP</b>  | Piperacillin                           |
| <b>P/T4</b> | Piperacillin /Tazobactam constant 4    |
| <b>TET</b>  | Tetracycline                           |
| <b>POS</b>  | Positive Control                       |

# One Health

**“Antimicrobial resistance is one of the greatest health challenges of our time, and we cannot leave it for our children to solve. Now is the time to forge new, cross-sector partnerships that will protect the medicines we have and revitalize the pipeline for new ones.”**

— Dr. Tedros Adhanom Ghebreyesus, WHO Director-General

The One Health approach is the One Health Quadripartite's global initiative to improve alignment across human, animal, and environmental sectors and, ultimately, tackle antimicrobial resistance.<sup>1</sup> Our unique combination of antimicrobial susceptibility testing solutions address human, animal, and environmental sectors. At Thermo Fisher Scientific we uphold the importance of this initiative as we understand that our health is all connected.

## Our health is all connected

Whether the patient has fur, feathers, or skin, accurate testing is paramount in diagnosing and treating infections.



- Our **clinical** portfolio helps optimize treatment for the critically ill, promoting better patient outcomes.
- Our **veterinary** solutions can be used to detect resistance in a wide range of animals, from pets and livestock to wildlife, using species-specific formats.
- Our **environmental** capabilities help to determine pathogen antimicrobial susceptibility patterns in food, animals, water and the environment, making us the system of choice for many global surveillance organizations.

**Ref:** 1. World Health Organization, One Health [website], <https://www.who.int/news-room/questions-and-answers/item/one-health>, accessed (21 December 2023).

The products are for Invitro Diagnostics use only.

For other details regarding the product, please refer to its IFU. INMBD1/04-25

©2025 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

The Thermo Fisher Scientific Logo is the sole intellectual property of Thermo Fisher Scientific Inc., USA

### Additional Information Available on Request:

Thermo Fisher Scientific India Pvt. Ltd.  
403-404, Delphi 'B' Wing, Hiranandani Business Park,  
Powai, Mumbai 400 076, India.  
customer.reach@thermofisher.com  
Call India toll-free: 1-800-222-230

Scan QR Code

